



PLEASE NOTE THAT THE MEASUREMENTS GIVEN IN THIS DRAWING ARE FOR YOUR INFORMATION ONLY

<b>POSS:</b> Qty. Description If the product is required to be manufactured differently, this must be approved by YONG LIFE-SAVING EQUIPMENT A/S and a specification of change must be filed in if the difference is approved. Unless otherwise specified the product must always be clearly marked with YONG part no., manufacturer ID and manufacturing year.		<b>Material:</b>	
<b>REVISION NO:</b> 00 <b>REVISION DESCRIPTION:</b>	<b>STATUS:</b>	<b>REVISION DATE:</b>	<b>APPROVAL INSTITUTE:</b>
<b>DATE:</b>	<b>DATE:</b>	<b>DATE:</b>	<b>DATE:</b>
<b>UNIT:</b>	<b>UNIT:</b>	<b>UNIT:</b>	<b>UNIT:</b>
<b>BASE MATERIAL:</b>	<b>SERVICE TREATMENT:</b>	<b>PRODUCTION TYPE:</b>	<b>DESCRIPTION:</b>
YONG LIFE-SAVING EQUIPMENT A/S Bøvlingsvej 13, DK-4170 Slønsbo V Tel.: +45 7611 8100 Fax: +45 7611 8101 <a href="mailto:info@yong-life.com">info@yong-life.com</a> <a href="http://www.yong-life.com">www.yong-life.com</a>	FIRST ANGLE PROJECTION GENERAL TOLERANCE:	Sales specification for container Part No.: 1033599 (016-9109) - W/POCKET 1033629 100DKS "A" and 150DKS "A" (Dossier) + "B"	DRAWING NO.: 42003064